




CONNECT


PROVIDER PORTAL

GET REPORTS | VIEW IMAGES | ORDER EXAMS | MAKE APPOINTMENTS

Quickly see completed exams from the previous 90 days.




Filter worklists in ways that make sense to you.




Exam Status

See exams by status—completed, scheduled, or in progress.



Filter by Patient Name

Enter the first or last name to receive a quick-pick list of matching results.




Filter by Exam Type

Click to filter exams by modality type.



Order Status

See orders by status—ordered or pending.



Filter by Ordering Provider

Filter ordering providers from the display by checking or unchecking them in the drop-down menu.



Order and Schedule

Order and schedule for most types of exams. Track orders by status. Access our appointment book and select available time-slots—in real time—including weekend and evening appointments.



Search

Search for your patients, or patients belonging to other providers in your practice, using common attributes. You can also search for patients outside your practice.



Get Help


If you need help getting started, solving technical issues, or learning how to use the portal, contact the portal's Web Team.



My Folder

Add exams to custom folders for convenient access in the future. Name your folders in ways that make sense to you to help you organize an exam archive.

View patient's exam history—providing access to reports and images.



Download Report
Download a PDF-copy of the final report.



Print Report
Print the report using the browser's print function.



Add to Folder
Add exams to custom folders for convenient access in the future.



Create Order
Create a new order for the selected patient.

The Images tab shows all images acquired during the patient's exam, organized in series.



Three Image Viewer Options

1. Quick Viewer

Click the lower-right corner of the thumbnail, indicated by a green highlight, to view a large, static display of the selected image.

2. Web Viewer

Access a basic set of controls to view and manipulate images with the Web Viewer, which works on any device.

3. Full Viewer

Use a full-featured viewer that displays diagnostic-quality images. Full Viewer must be installed on your computer prior to first use.



Account

The Account Page allows you to define settings that are specific to your portal account:

My Profile

- Update your personal information, including email and phone numbers.

Settings

- Refine your view preferences for how patients are displayed on the Home Page and Orders Page.

Security

- Change your password.
- Change your security answers. These are helpful to retrieve your user ID or password, if you forget in the future.

Orders

- Change the default settings.
- Authorized others to order and schedule on your behalf (ordering providers only).
- See which providers you can order and schedule for without their review and approval.

Image Viewer

- Define the default Image Viewer.